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Aptal – Resin Root Filling

Composition:

powder contains: zinc oxide, silica (aerosil), rosin, silver powder.
The liquid contains eugenol, rosin, olive oil, turpentine oil and Peruvian balsam.

Indications:

To fill root canals in teeth following root canal treatment.

Contraindications:

No systemic side effects are known to date.

Side effects/Interactions:

In individual cases, allergic reactions have been described. Direct contact with synthetic filling materials (composites) must be avoided.

Dosage instructions, mode and duration of application:

The root canal filling material is introduced into the root canal using appropriate instruments. The root canal must be dried beforehand because the presence of moisture makes the material set more quickly. The use of posts is common scientific knowledge in endodontics. The powder is gradually wetted with the mixing fluid on the smooth side of as thick a glass plate as possible to facilitate even distribution of the paste introduced into the root canal using a spiral instrument (lentula). One of the principal advantages is that the paste adheres very firmly to the canal walls.

Once the paste becomes viscous, it is no longer advisable to work with the lentula because resistance will be too great and the spiral may fracture. The hardening time of the Aptal resin root filling does not depend so much on the mixing ratio; the temperature has a much greater effect. To slow down setting at higher ambient temperatures, chilled glass plates of the type also recommended for use when working with silicate cements are used. It is pointless to mix the Aptal resin root filling until it is very thin in an attempt to achieve slower setting. The normal thick, creamy consistency ensures as tightly packed a root canal filling as possible. Gutta-percha points coated with paste are then introduced. The coronal region of the root canal filling is covered with the rest of the increasingly viscous paste, which can, if necessary, be thickened by the addition of more powder. This protective layer, which is more reliable than a phosphate cement underlay, is pressed down using a cotton-wool pellet slightly moistened with alcohol. If a stronger filling is required, a suitable pear-shaped or spherical filling instrument is used, rolled beforehand in the powder residue on the glass plate or greased lightly with Vaseline, to prevent the viscous paste from adhering to the instrument. Alternatively, phosphate cements can be used, e.g. if a permanent synthetic filling is to be used. If the root canal is overfilled and/or if the root filling material enters the mandibular canal (Canalis mandibulae) or the maxillary sinus, the question of revision must be clarified.

Special instructions:

This medical product must not be used after the expiry date.
Keep out of reach of children.

Available in the following packs: Aptal resin root filling powder, Contents: 40 g Article No. 1003
Aptal resin root filling liquid, Contents: 10 ml Article No. 1004

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