SPEIKO® Carboxylzement (vormals Zhanelka)

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Stand der Information: Februar 2006

SPEIKO® Carboxylate Cement (previously Zhanelka)

Composition:

Powder contains : Zinc oxide, Silicon oxide, Aluminium oxide, Calcium oxide, Sodium oxide, Magnesium oxide

Liquid contains : Polyacrylic acid

Application:

SPEIKO® Carboxylate Cement is used as lining for fillings of amalgam, silicate, or synthetic material and for fixing inlays, crowns, and bridges.

Contraindications:

Because of the possibility of pulp injury the product may not be applied on the dentine close to the pulp or the open pulp. Using SPEIKO® Carboxylate cement for temporary filling, you will have ablated material as long as it is incorporated.

Abnormal sensitivity against one of the components.

Side effects / Interactions:

After using for lining or temporary filling pulp irritations may occur.

Dosage instructions, form and duration of application:

Mix powder with SPEIKO® Carboxylate Cement liquid in volume ratio 2,2 : 1, taking care that the powder is mixed into the liquid and that the mixing process does not exceed one minute. Then a suitable consistency should be obtained.

A glass plate or mixing pad is used for mixing. On the plate SPEIKO® Carboxylate Cement remains longer plastic when it is cooled.

The surfaces or teeth and cavity should be dry and clean, so cement may be filled immediately after mixing. The setting of the cement is completed within 2,5 - 8 min.

A right setting may be tested with the rest remaining on the mixing pad.

To achieve a good adhesion of SPEIKO® Carboxylate cement the teeth surfaces and restorations must be free of oil, grease, wax and other contaminations.

Special instructions:

The liquid contains polyacrylic acid and has a corrosive effect. In the event of inadvertent contact with the skin, mucosa or eyes, rinse immediately with plenty of water. Lack of appropriate pulp protection may result in pulp necrosis. The patient should not eat for the first 2 hours after receiving the filling. This medical product should not be used after the expiry date.

Keep out of reach of children.

Available in the following packs:

SPEIKO® Carboxylate Cement Powder	Contents: 100 g	Art.-no. 1104
SPEIKO® Carboxylate Cement liquid	Contents: 20 ml	Art.-no. 1105

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SPEIKO® - Dr. Speier GmbH P.O. Box 7204 D-48039 Münster