SPEIKOFIL-AQUA

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Artikel-Nr. 1470

Stand Januar 2000

SPEIKOFIL-AQUEOUS
glass ionomer cement mixable with water

DIN EN 29917

Composition:

Base materials: silicone dioxide, aluminium fluoride, iron oxide (colouring), polyacrylic acid.

Uses:

SPEIKOFIL aqueous is suitable as a base, for the treatment of primary teeth and as a build-up filler. It is indicated for tooth neck cavities (Black-class V), Black class III restorations (approximal anterior cavities) and also for the treatment of root caries.

Contraindications:

Because of possible damage to the pulp the product should not be applied to dentine near the pulp or to the open pulp. Use of calcium hydroxide as a base (Speikocal, Calasept). Hypersensitivity to the components.

Side effects and interaction with other materials:

Used in accordance with the directions for use there are at present no known side effects and no known interaction with other materials.

Measuring instructions, method of use and working time:

The mixing ratio is 1 measuring spoon to 1 drop of water (1:1). Best results are obtained with 5.12 gr. to 5.53 gr. powder to 1 gr. of water at 23°C and 50% humidity.

The powder should be well shaken before use.

The bottle should be half filled with de-ionized, distilled or tap water to facilitate accurate measuring .

Using level spoonfuls of the measuring spoon put the required amount of powder onto a glass slab (or mixing block) and divide into three portions. Drop on water and mix the first portion with a plastic spatula, then mix in the second and third portions. The mixing procedure should not take longer than 30 sec..

The working time after mixing is 2 min. Using the optimum mixing ratio setting begins after 4 to 5 min..

Ensure that all surfaces of the tooth and the cavity for the filling are prepared, cleaned and dried so that the cement can be applied immediately after mixing. The filling can be polished after ca. 7 min. .

If required it can be checked that the cement is setting properly using the residue remaining on the glass slab.

Caution:

Liquid contains polyacrylic acid. Avoid contact with the eyes.

Powder is hygroscopic. Close container tightly after use.

Without a suitable base pulp necrosis may occur.

The patient should not eat anything for 2 hours after insertion of the filling.

This pharmaceutical product should not be used after the expiry date.

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